

**Cat No:** Kab21891

**Product name:** Chikungunya E2 monoclonal antibody, clone KB583

**Pack Size:** 1 mg

**Host:** Mouse

**Isotype:** IgG1, kappa

**Purity:** >90%

**Storage:** Store at 4°C. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing.

**Specificity:** Viruses

**Type:** Unconjugated

**Immunogen:** Recombinant Chikungunya virus envelope 2 protein.

**Clone :** KB583

**Form:** Liquid

**Buffer:** In 15 mM Potassium Phosphate Buffer, pH 7.2 with 0.85% NaCl (0.05% sodium azide).

**Applications:** ELISA, LFIA

*Note: Optimal concentrations/dilutions should be determined by the end user.*

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

OUR KINESISDX COMMITMENT: GUARANTEED QUALITY WITH EXPERT TECHNICAL SUPPORT