

Cat No: Kab21892

Product name: Chikungunya E2 monoclonal antibody, clone KB584

Pack Size: 1 mg

Host: Mouse

Isotype: IgG2b, kappa

Purity: >90%

Storage: Store at 4°C. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing.

Specificity: Mouse

Type: Unconjugated

Immunogen: Recombinant Chikungunya virus envelope 2 protein.

Clone : KB584

Form: Liquid

Buffer: In 15 mM Potassium Phosphate Buffer, pH 7.2 with 0.85% NaCl (0.05% sodium azide).

Applications: ELISA, LFIA

Note: Optimal concentrations/dilutions should be determined by the end user.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

OUR KINESISDX COMMITMENT: GUARANTEED QUALITY WITH EXPERT TECHNICAL SUPPORT