

Cat No: KDXC529037

Product Name: Rabies Virus Glycoprotein Mouse anti-Rabies Virus Monoclonal (Biotin) Antibody

Pack Size: 100 µL

Host: Mouse

Isotype: IgG2a

Storage: Short term: store at 4°C. Long term: store at -20°C. Avoid repeated freeze-thaw cycles.

Description: Rabies Virus Glycoprotein antibody is a biotin-conjugated mouse monoclonal antibody to rabies virus Glycoprotein.

Purification: Protein A purified

Specificity: Recognizes glycoprotein of rabies virus

Reactivity: Rabies Virus (tested or 100% immunogen sequence identity)

Immunogen: Purified rabies virus

Tested Applications:

ELISA

Formulation: PBS, pH 7.2, No preservatives added, No stabilizing proteins added

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

OUR KINESISDX COMMITMENT: GUARANTEED QUALITY WITH EXPERT TECHNICAL SUPPORT