

**Cat No: KDXC529038**

**Product Name:** Rabies Virus Glycoprotein Mouse anti-Rabies Virus Monoclonal (Biotin) Antibody

**Pack Size:** 100 µL

**Host:** Mouse

**Isotype:** IgG2a

**Storage:** Short term: store at 4°C. Long term: store at -20°C. Avoid repeated freeze-thaw cycles.

**Description:** Rabies Virus Glycoprotein antibody is a biotin-conjugated mouse monoclonal antibody to rabies virus Glycoprotein.

**Purification:** Protein A purified

**Specificity:** Recognizes glycoprotein of rabies virus

**Reactivity:** Rabies Virus (tested or 100% immunogen sequence identity)

**Immunogen:** Purified rabies virus

**Tested applications:**

ELISA

**Formulation:** PBS, pH 7.2, No preservatives added, No stabilizing proteins added

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

OUR KINESISDX COMMITMENT: GUARANTEED QUALITY WITH EXPERT TECHNICAL SUPPORT